



Orthopedic implants



Catalogue
New Collection




LOCKING PLATES	1
LOCKING SCREW	16
BONE PLATE	18
BONE SCREW	30
INTERLOCKING NAIL SYS	37
INTERLOCKING BOLT SYS	42
PROSTHESIS SYSTEM	44
SPINE SYSTEM	46



SMALL COMPRESSION LOCKING PLATE 3.5MM

Catalog no.	011.155.03 to 011.155.16
Screw \varnothing	3.5mm
Material	Stainless Steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



LCP SMALL T PLATE (OBLIQUE) 3.5MM - LEFT AND RIGHT

Catalog no.	011.151.02 to 011.151.10
Screw \varnothing	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	2,3,4,5,6,7,8,9,10
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



LCP SMALL T PLATE (RIGHT ANGLE) 3.5MM

Catalog no.	011.150.02 to 011.150.10
Screw \varnothing	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	2,3,4,5,6,7,8,9,10
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)

LCP PROXIMAL HUMERUS PLATE 3.5MM (TYPE 1)

Catalog no.	011.156.04 to 011.156.12
Screw ø	3.5mm & 4.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4,6,8,10,12
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



LCP PROXIMAL HUMERUS PLATE 3.5MM (TYPE 2)

Catalog no.	011.1910.03 to 011.1910.10
Screw ø	3.5mm & 4.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



LCP PROXIMAL HUMERUS PLATE 3.5MM (TYPE 3)

Catalog no.	011.1911.03 to 011.1911.08
Screw ø	3.5mm & 4.0mm
Material	Stainless steel 316L, Titanium
Holes	3,4,5,6,7,8
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)





LCP CLOVERLEAF PLATE 3.5MM

Catalog no.	011.154.03 to 011.154.03
Screw \varnothing	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13
Proximal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



LCP RECONSTRUCTION PLATE 3.5MM ROUND HOLE

Catalog no.	011.153.03 to 011.153.18
Screw \varnothing	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)




LCP RECONSTRUCTION PLATE 5.0MM (ROUND HOLE)

Catalog no.	011.1954.03 to 011.1954.18
Screw \varnothing	5.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (Diff. 2.0mm)

LCP RECONSTRUCTION PLATE 3.5MM (DYNAMIC HOLE)

Catalog no.	011.1912.03 to 011.1912.16
Screw ø	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16
Shaft Part	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



LCP RECONSTRUCTION PLATE 5.0MM (DYNAMIC HOLE)

Catalog no.	011.1913.03 to 011.1913.16
Screw ø	5.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (Diff. 2.0mm)



DISTAL TIBIA LOCKING PLATE LEFT & RIGHT

Catalog no.	011.157.04 to 011.157.12
Screw ø	5.0mm & 4.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4, 6, 8, 10, 12,
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking screw 5.0mm
Screw Length	
(Locked Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)





MEDIAL DISTAL TIBIA LOCKING PLATE 3.5MM LEFT AND RIGHT

Catalog no.	011.1909.05 to 011.1909.13
Screw ø	3.5mm & 4.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	5, 7, 9, 11, 13
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



DISTAL TIBIA ANATOMICAL LOCKING PLATE - LEFT & RIGHT

Catalog no.	011.1001.04 to 011.1001.14
Screw ø	3.5mm & 4.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4, 6, 8, 10, 12, 14
Proximal Screw	Locking Screw 4.0mm
Screw Length	20mm to 70mm (Diff. 5.0mm)
Distal Screw	Locking Screw 3.5mm
Screw Length	10mm to 40mm (Diff. 2.0mm)



NARROW LOCKING COMPRESSION PLATE 5.0MM

Catalog no.	011.166.04 to 011.166.18
Screw ø	5.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Shaft Part	Locking Screw 5.0mm
Screw Length	
(Locked Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)

BROAD LOCKING COMPRESSION PLATE 4.8MM

Catalog no. 011.167.04 to 011.167.20
Screw ø 5.0mm
Material Stainless steel 316L,316LVM, Titanium
Holes 4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19,20
Shaft Part Locking Screw 5.0mm
Screw Length
(Locked Screw 5.0mm) 10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5.0mm) 55mm to 90mm (Diff. 5.0mm)



L BUTTRESS LOCKING PLATE LEFT & RIGHT

Catalog no. 011.148.04 to 011.148.12
Screw ø 5.0mm
Material Stainless steel 316L,316LVM, Titanium
Holes 4,5,6,7,8,9,10,11,12
Proximal Screw
Screw Length
Distal Screw Locking Screw 5.0mm
Screw Length
(Locked Screw 5.0mm) 10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5.0mm) 55mm to 90mm (Diff. 5.0mm)



T PLATE LOCKING 5.0MM

Catalog no. 011.162.04 to 011.162.12
Screw ø 5.0mm
Material Stainless steel 316L,316LVM, Titanium
Holes 4,5,6,7,8,9,10,11,12
Proximal Screw
Screw Length
Distal Screw Locking Screw 5.0mm
Screw Length
(Locked Screw 5.0mm) 10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5.0mm) 55mm to 90mm (Diff. 5.0mm)





T BUTTRESS LOCKING PLATE 5.0MM

Catalog no.	011.1262.04 to 011.1262.12
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12
Proximal Screw	
Screw Length	
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locked Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)



DISTAL HUMERUS LOCKING PLATE LEFT & RIGHT

Catalog no.	011.1925.03 to 011.1925.15
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,15
Shaft Part	Locking Screw 3.5mm
Screw Length	10.0mm to 60.0mm (diff. 2.0mm)



DISTAL HUMERUS WITH LATERAL SUPPORT LOCKING PLATE LEFT & RIGHT

Catalog no.	011.1907.03 to 011.1907.15
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,15
Shaft Part	Locking Screw 3.5mm
Screw Length	10.0mm to 60.0mm (diff. 2.0mm)

EXTRA-ARTICULAR DISTAL HUMERUS LOCKING PLATE LEFT & RIGHT

Catalog no.	011.1941.04 to 011.1941.14
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,6,8,10,12,14
Shaft Part	Locking Screw 3.5mm
Screw Length	10.0mm to 60.0mm (diff. 2.0mm)



LATERAL TIBIAL LOCKING PLATE 5.0MM/6.5MM – LEFT & RIGHT

Catalog no.	011.161.03 to 011.161.12
Screw ø	5.0mm & 6.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12
Proximal Screw	Locking Screw 6.5mm
Screw Length	30mm to 95mm (Diff. 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locking Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locking Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)



LATERAL TIBIAL LOCKING PLATE 5.0MM – LEFT & RIGHT

Catalog no.	011.159.03 to 011.159.13
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13
Proximal Screw	Cancellous Locking Screw 5.0mm
Screw Length	30mm to 80mm (Diff. 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locking Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locking Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)





DISTAL FEMUR LOCKING PLATE 5.0MM/6.5MM (LEFT AND RIGHT)

Catalog no.	011.163.03 to 011.163.13
Screw ø	5.0mm & 6.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13
Proximal Screw	Cancellous Locking Screw 6.5mm
Screw Length	30mm to 95mm (Diff. 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locked Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locked Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)



DISTAL FEMUR LOCKING PLATE LEFT & RIGHT

Catalog no.	011.164.03 to 011.164.13
Screw ø	5.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13
Proximal Screw	Cancellous Locking Screw 5.0mm
Screw Length	30mm to 80mm (Diff. 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locking Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locking Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)



CONDYLAR BUTTERS LOCKING PLATE LEFT & RIGHT

Catalog no.	011.176.03 to 011.176.13
Screw ø	5.0mm & 6.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13
Proximal Screw	Cancellous Locking Screw 6.5mm
Screw Length	30mm to 95mm (Diff. 5.0mm)
Distal Screw	Locking Screw 5.0mm
Screw Length	
(Locking Screw 5.0mm)	10mm to 50mm (Diff. 2.0mm)
(Locking Screw 5.0mm)	55mm to 90mm (Diff. 5.0mm)

DHS PLATE LOCKING LONG BARRAL 130° /135°

Catalog no.	011.165.02 to 011.165.12
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	2,3,4,5,6,7,8,9,10,11,12
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (Diff. 2.0mm) 55mm to 90mm (Diff. 5.0mm)



DHS PLATE LOCKING SHORT BARRAL 130° /135°

Catalog no.	011.152.02 to 011.152.12
Screw ø	5.0mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	2,3,4,5,6,7,8,9,10,11,12
Shaft Part	Locking Screw 5.0mm
Screw Length	10mm to 50mm (Diff. 2.0mm) 55mm to 90mm (Diff. 5.0mm)



D.C.S PLATE LOCKING 90°

Catalog no.	011.1937.03 to 011.1937.16
Screw ø	5.0mm & 6.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12,13,14,15,16
Locking Screw 5.0mm	10mm to 50mm (Diff. 2.0mm) 55mm to 90mm (Diff. 5.0mm)
Cancellous Locking	
Screw 6.5mm	30mm to 95mm (Diff. 5.0mm)





DISTAL FEMUR LOCKING PLATE 5.0MM/6.5MM (LEFT AND RIGHT)

Catalog no.	011.1928.04 to 011.1928.08
Screw ø	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	04, 06, 08
Locking Screw	
3.5mm	12mm to 40mm (Diff. 2.0mm)



PROXIMAL FEMUR LOCKING PLATE

Catalog no.	011.1938.05 to 011.1938.13
Screw ø	5.0mm & 6.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	05, 07, 09, 11, 13
Locking Screw	
5.0mm	10mm to 50mm (Diff. 2.0mm) 55mm to 90mm (Diff. 5.0mm)
Cancellous	
Screw 6.5mm	30mm to 100mm (Diff. 5.0mm)



PROXIMAL TIBIA LOCKING PLATE LEFT & RIGHT

Catalog no.	011.1986.04 to 011.1986.14
Screw ø	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	04, 06, 08, 10, 12, 14
Locking Screw	
3.5mm	12mm to 40mm (Diff. 2.0mm)


LCP TOMIFIX MID HIGH T PLATE

Catalog no.	011.1252.03 & 011.1252.04
Screw ø	5.0mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	03, 04
Locking Screw	
5.0mm	18mm to 50mm (Diff. 2.0mm) 55mm to 90mm (Diff. 5.0mm)



CALCANEUS LOCKING PLATE LEFT & RIGHT

Catalog no.	011.1940.70mm & 011.1940.75mm
Screw ø	3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	70mm, 75mm
Locking Screw	
3.5mm	12mm to 40mm (Diff. 2.0mm)



DISTAL TIBIA WITH TAB LOCKING PLATE - LEFT & RIGHT

Catalog no.	011.1967.05 to 011.1967.13
Screw ø	4.0 mm & 3.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	05, 07, 09, 11, 13
Locking Screw	
3.5mm	10mm to 40mm (Diff. 2.0mm)
Locking Screw	
4.0mm	20mm to 70mm (Diff. 5.0mm)





2.4MM DISTAL RADIUS LOCKING PLATE STRAIGHT

Catalog no.	011.193.05 to 011.193.13
Screw ø	2.4mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	5, 7, 9, 11, 13
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)



2.4MM L DISTAL RADIUS LOCKING PLATE (HEAD 2HOLE) – LEFT & RIGHT

Catalog no.	011.191.3 & 011.191.4
Screw ø	2.4mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	2+3, 2+4
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)



2.4MM L DISTAL RADIUS LOCKING PLATE (HEAD 3HOLE) – LEFT & RIGHT

Catalog no.	011.190.3 & 011.190.4
Screw ø	2.4mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3+3, 3+4
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)

2.4MM L DISTAL RADIUS OBLIQUE LOCKING PLATE – LEFT & RIGHT

Catalog no. 011.194. 3 & 011.194. 4
Screw ø 2.4mm
Material Stainless steel 316L,316LVM, Titanium
Holes 3+3, 3+4
Locking Screw
2.4mm 06mm to 30mm (Diff. 2.0mm)



2.4MM T DISTAL RADIUS LOCKING PLATE (HEAD 3HOLE)

Catalog no. 011.192. 3 & to 011.192. 4
Screw ø 2.4mm
Material Stainless steel 316L, 316LVM, Titanium
Holes 3+3, 3+4
Locking Screw
2.4mm 06mm to 30mm (Diff. 2.0mm)



2.4MM DISTAL RADIUS LOCKING PLATE (HEAD 5HOLE) – LEFT & RIGHT

Catalog no. 011.197. 3 & to 011.197. 5
Screw ø 2.4mm
Material Stainless steel 316L,316LVM, Titanium
Holes 5+3, 5+5
Locking Screw
2.4mm 06mm to 30mm (Diff. 2.0mm)





2.4 MM BUTTRESS LOCKING PLATE (HEAD 5HOLE) – LEFT & RIGHT

Catalog no.	011.198
Screw \varnothing	2.4mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	5+3
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)



TWO-COLUMN DISTAL RADIUS PLATE LEFT & RIGHT

Catalog no.	011.1946.02 to 011.1946.04
Screw \varnothing	2.4mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	02, 03, 04
Locking Screw	
2.4mm	06mm to 30mm (Diff. 2.0mm)

LCP SCREW DIA 2.4 MM

Catalog no	011.286.06 to 011.286.30
Thread Diameter	2.4mm
Use drills Diameter	2.0mm
Drive	hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	2.7mm
Length	06mm to 30mm (Diff. 2.0mm)




LCP SCREW DIA 2.7 MM

Catalog no	011.204.06 to 011.204.30
Thread Diameter	2.7mm
Use drills Diameter	2.0mm
Drive	hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	2.7mm
Length	06mm to 30mm (Diff. 2.0mm)



LCP SCREW DIA 3.5 MM

Catalog no	011.268.10 to 011.268.50
Thread Diameter	
Use drills Diameter	2.7mm
Drive	Hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	3.5mm
Length	10mm to 50mm (Diff. 2.0mm)



LCP SCREW DIA 4.0 MM

Catalog no	011.270.25 to 011.270.70
Thread Diameter	4.0mm
Use drills Diameter	3.2mm
Drive	Hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	4.0mm
Length	25mm to 70mm (Diff. 2.0mm)




LOCKING SCREW



LCP SCREW DIA 5.0 MM

Catalog no	011.269.16 to 011.269.95
Thread Diameter	5.0mm
Use drills Diameter	4.0mm
Drive	Hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	5.0mm
Length	16mm to 50mm (Diff. 2.0mm) 55mm to 95mm (Diff. 5.0mm)



LCP SCREW CANCELLOUS DIA 5.0 MM

Catalog no	011.2930.30 to 011.2930.80
Thread Diameter	5.0mm
Use drills Diameter	4.0mm
Drive	Hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	5.0mm
Length	30mm to 50mm (Diff. 2.0mm) 55mm to 80mm (Diff. 5.0mm)



LCP SCREW DIA 6.5 MM

Catalog no	011.271.35 to 011.271.95
Thread Diameter	6.5mm
Use drills Diameter	4.5mm
Drive	hex
Material	Stainless steel 316L, 316LVM, Titanium
Size	6.5mm
Length	35mm to 95mm (Diff. 5.0mm)

RECONSTRUCTION PLATE FOR 2.7MM SCREW ROUND HOLE

Catalog no.	011.300.04 to 011.300.14
Screw ø	2.7mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,14
Screw	Cortical Screw 2.7mm
Screw Length	06mm to 30mm (Diff. 2.0mm)



DYNAMIC COMPRESSION PLATE FOR 2.7 MM SCREW

Catalog no.	011.320.04 to 011.320.14
Screw ø	2.7mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,14
Screw	Cortical Screw 2.7mm
Screw Length	06mm to 30mm (Diff. 2.0mm)



ONE THIRD TUBULAR PLATE FOR 3.5 MM SCREW

Catalog no.	011.321.04 to 011.321.14
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,14
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)





RECONSTRUCTION PLATE FOR 3.5 MM SCREW ROUND HOLE

Catalog no.	011.322.04 to 011.322.18
Screw \varnothing	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,13,14,15,16,17,18
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



RECONSTRUCTION PLATE FOR 3.5 MM SCREW DC HOLE

Catalog no.	011.3221.04 to 011.3221.14
Screw \varnothing	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,14
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



RECONSTRUCTION CURVED PLATE FOR 3.5MM SCREW

Catalog no.	011.301.04 to 011.301.16
Screw \varnothing	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,08,09,10,11,12,13,14,15,16
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)

SMALL DCP PLATE FOR 3.5 MM SCREW

Catalog no.	011.325.04 to 011.325.14
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,14
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



LC DCP PLATE SMALL FOR 3.5 MM SCREW

Catalog no.	011.333.04 to 011.333.14
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06,07,08,09,10,11,12,14
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



SMALL T PLATE FOR 3.5 MM SCREW

Catalog no.	011.326.03 to 011.326.10
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	03,04,05,06,07,08,09,10
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)





SMALL T OBLIQUE PLATE FOR 3.5 MM SCREW

Catalog no.	011.327.03 to 011.327.10
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	03,04,05,06,07,08,09,10
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



ACRO PLATE FOR 3.5 MM SCREW LEFT & RIGHT

Catalog no.	011.3931.04 to 011.3931.06
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	04,05,06
Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)



RECONSTRUCTION PLATE FOR 4.5 MM SCREW ROUND HOLE

Catalog no.	011.323.04 to 011.323.18
Screw ø	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)

RECONSTRUCTION PLATE FOR 4.5 MM SCREW DC HOLE

Catalog no.	011.324.04 to 011.324.14
Screw \varnothing	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14,
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



SEMI TUBULAR PLATE FOR 4.5 MM SCREW

Catalog no.	011.3225.04 to 011.3225.14
Screw \varnothing	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14,
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



NARROW DCP PLATE FOR 4.5 MM SCREW

Catalog no.	011.337.04 to 011.337.14
Screw \varnothing	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)





LC DCP PLATE NARROW FOR 4.5 MM SCREW

Catalog no.	011.334.04 to 011.334.14
Screw \varnothing	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



BROAD DCP PLATE FOR 4.5 MM SCREW

Catalog no.	011.338.04 to 011.338.18
Screw \varnothing	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



LC DCP PLATE BROAD FOR 4.5 MM SCREW

Catalog no.	011.358.06 to 011.358.14
Screw \varnothing	4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	6,7,8,9,10,11,12,14
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)

CLOVERLEAF PLATE

Catalog no. 011.328.03 to 011.328.12
Screw ø
Material Stainless steel 316L,316LVM, Titanium
Holes 3,4,5,6,7,8,9,10,11,12
Screw
Screw Length



T PLATE FOR 4.5 MM SCREW

Catalog no. 011.329.03 to 011.329.12
Screw ø 6.5mm & 4.5mm
Material Stainless steel 316L,316LVM, Titanium
Holes 3,4,5,6,7,8,9,10,11,12
Proximal Screw Cancellous Screw 6.5mm
Screw Length 30mm to 110mm (Diff. 5.0mm)
Distal Screw Cortical Screw 4.5mm
Screw Length 16mm to 70mm (Diff. 2.0mm)



T BUTTRESS PLATE FOR 4.5 MM SCREW

Catalog no. 011.330.03 to 011.330.12
Screw ø 6.5mm & 4.5mm
Material Stainless steel 316L,316LVM, Titanium
Holes 3,4,5,6,7,8,9,10,11,12
Proximal Screw Cancellous Screw 6.5mm
Screw Length 30mm to 110mm (Diff. 5.0mm)
Distal Screw Cortical Screw 4.5mm
Screw Length 16mm to 70mm (Diff. 2.0mm)





L BUTTRESS PLATE FOR 4.5 MM SCREW LEFT & RIGHT

Catalog no.	011.332.03 to 011.332.12
Screw ø	6.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	3,4,5,6,7,8,9,10,11,12
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



LATERAL TIBIAL PLATE (HOCKEY) FOR 4.5 MM SCREW LEFT & RIGHT

Catalog no.	011.335.04 to 011.335.14
Screw ø	6.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



CONDYLAR BUTTRESS PLATE FOR 4.5 MM SCREW LEFT & RIGHT

Catalog no.	011.336.04 to 011.336.14
Screw ø	6.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)

DHS PLATE (DCP HOLES) SHORT BARREL

Catalog no.	011.339.04 to 011.339.14
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14
Degree	120°, 125°, 130°, 135°, 140°, 145°
Proximal Screw	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



DHS PLATE (ROUND HOLES) SHORT BARREL

Catalog no.	011.339.04R to 011.339.18R
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Degree	120°, 125°, 130°, 135°, 140°, 145°
Proximal Screw	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



DHS PLATE (DCP HOLES) LONG BARREL

Catalog no.	011.340.04 to 011.340.14
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,14
Degree	120°, 125°, 130°, 135°, 140°, 145°
Proximal Screw	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)





DHS PLATE (ROUND HOLES) LONG BARREL

Catalog no.	011.340.02R to 011.340.18R
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18
Degree	120°, 125°, 130°, 135°, 140°, 145°
Proximal Screw	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



DHS PLATE (ROUND HOLES) PEDIATRICS

Catalog no.	011.331.04 to 011.331.10
Screw ø	
Material	Stainless steel 316L, 316LVM, Titanium
Holes	4,5,6,7,8,9,10
Degree	125°, 130°, 135°
Proximal Screw	Lag Screw
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)



95° DCS PLATE (DCP HOLES) TOTAL HOLE

Catalog no.	011.341.06 to 011.341.20
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L, 316LVM, Titanium
Holes	6,7,8,9,10,11,12,13,14,15,16,17,18,19,20
Proximal Screw	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)

DCS PLATE (ROUND HOLES) 95°

Catalog no.	011.342.06 to 011.342.20
Screw ø	12.5mm & 4.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	6,7,8,9,10,11,12,13,14,15,16,17,18,19,20
Proximal Screw	Lag Screw 12.5mm
Screw Length	50mm to 115mm (Diff. 5.0mm)
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)



Y RECONSTRUCTION PLATE FOR 3.5mm SCREW

Catalog no.	011.302.01 to 011.302.10
Screw ø	3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	1,2,3,4,5,6,7,8,9,10
Distal Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 4.0mm
Screw Length	10mm to 30mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 4.0mm
Screw Length	35mm to 70mm (Diff. 5.0mm)



CALCANEAL PLATE 3.5MM

Catalog no.	011.310
Screw ø	4.0mm & 3.5mm
Material	Stainless steel 316L,316LVM, Titanium
Holes	Small , Medium , Large
Distal Screw	Cortical Screw 3.5mm
Screw Length	10mm to 60mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 4.0mm
Screw Length	10mm to 30mm (Diff. 2.0mm)
Distal Screw	Cancellous Screw 4.0mm



HOOK PLATE FOR 3.5mm SCREW

Catalog no.	011.3944.04 to 10
Screw \varnothing	3.5mm
Material	Stainless steel 316L,316LVM ,Titanium
Holes	3,4,5,6,7,8,9,10,11,12
Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)
Screw	Cortical Screw 4.5mm



COBRA HEAD PLATE FOR 4.5mm SCREW LEFT & RIGHT

Catalog no.	011.311.04 to 011.311.14
Screw \varnothing	
Material	Stainless steel 316L,316LVM, Titanium
Holes	4,5,6,7,8,9,10,11,12,13,14
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)

TROCENTRIC PLATE FOR DHS

Catalog no.	011.3948
Screw \varnothing	
Material	Stainless steel 316L,316LVM, Titanium
Holes	4 Hole
Distal Screw	Cortical Screw 4.5mm
Screw Length	16mm to 70mm (Diff. 2.0mm)
Proximal Screw	Cancellous Screw 6.5mm
Screw Length	30mm to 110mm (Diff. 5.0mm)

EIGHT PLATE

Catalog no.	011.366
Screw \varnothing	
Material	Stainless steel 316L,316LVM, Titanium
Holes	4 Hole
Distal Screw	Cortical Screw 4.5mm



CORTICAL SCREW (SELF TAPPING) 1.5MM

Catalog no.	011.49.06 to 011.49.30
Material	Stainless steel 316L,316LVM, Titanium
Size	1.5mm
Length	06mm to 30mm (Diff. 2.0mm)



CORTICAL SCREW 2.0MM

Catalog no.	011.49.06 to 011.49.30
Material	Stainless steel 316L,316LVM, Titanium
Size	1.5mm
Length	06mm to 30mm (Diff. 2.0mm)



CORTICAL SCREW 2.4MM

Catalog no.	011.485.06 to 011.485.30
Material	Stainless steel 316L,316LVM, Titanium
Size	2.4mm
Length	06mm to 30mm (Diff. 2.0mm)



CORTICAL SCREW 2.7MM

Catalog no.	011.401.06 to 011.401.30
Material	Stainless steel 316L,316LVM, Titanium
Size	2.7mm
Length	06mm to 30mm (Diff. 2.0mm)





CORTICAL SCREW (14 TPI) 3.5MM

Catalog no.	011.402.10 to 011.402.30
Material	Stainless steel 316L,316LVM, Titanium
Size	3.5mm
Length	10mm to 60mm (Diff. 2.0mm)



CORTICAL SCREW (14 TPI) SELF TAPPING 3.5MM

Catalog no.	011.01.402.10 to 011.01.402.60
Material	Stainless steel 316L,316LVM, Titanium
Size	3.5mm
Length	10mm to 60mm (Diff. 2.0mm)



CORTICAL SCREW 3.5MM (20 TPI)

Catalog no.	011.403.10 to 011.403.10.60
Material	Stainless steel 316L,316LVM Titanium
Size	3.5mm
Length	10mm to 60mm (Diff. 2.0mm)



CORTICAL SCREW (20 TPI) SELF TAPPING 3.5MM

Catalog no.	011.01.403.10 to 011.01.403.60
Material	Stainless steel 316L,316LVM, Titanium
Size	3.5mm
Length	10mm to 60mm (Diff. 2.0mm)

CANCELLOUS SCREW (SHORT THREAD) 4.0MM

Catalog no. 011.404.10 to 011.404.70
Material Stainless steel 316L,316LVM, Titanium
Size 4.0mm
Length 10mm to 30mm(Diff. 2.0mm)
 35mm to 70mm(Diff. 5.0mm)



CANCELLOUS SCREW (FULL THREAD) 4.0MM

Catalog no. 011.405.10 to 011.405.70
Material Stainless steel 316L,316LVM, Titanium
Size 4.0mm
Length 10mm to 30mm(Diff. 2.0mm)
 35mm to 70mm(Diff. 5.0mm)



CANNULATED SCREW (SHORT THREAD) 4.0MM

Catalog no. 011.406.10 to 011.406.70
Material Stainless steel 316L,316LVM, Titanium
Size 4.0mm
Length 10mm to 30mm(Diff. 2.0mm)
 35mm to 70mm(Diff. 5.0mm)



CANNULATED SCREW (FULL THREAD) 4.0MM

Catalog no. 011.407.10 to 011.407.70
Material Stainless steel 316L,316LVM, Titanium
Size 4.0mm
Length 10mm to 30mm(Diff. 2.0mm)
 35mm to 70mm(Diff. 5.0mm)





CORTICAL SCREW 4.5MM DIA

Catalog no.	011.409.16 to 011.409.70
Material	Stainless steel 316L,316LVM, Titanium
Size	4.5mm
Length	16mm to 70mm(Diff. 2.0mm)



CORTICAL SCREW (SELF TAPPING) 4.5MM

Catalog no.	011.01.409.12 to 011.01.409.70
Material	Stainless steel 316L,316LVM, Titanium
Size	4.5mm
Length	16mm to 70mm(Diff. 2.0mm)



MALLEOLAR SCREW 4.5MM

Catalog no.	011.410.20 to 011.410.100
Material	Stainless steel 316L,316LVM, Titanium
Size	4.5mm
Length	20mm to 100mm(Diff. 5.0mm)



CANCELLOUS SCREW (16TH) 6.5MM

Catalog no.	011.411.20 to 011.411.110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	20mm to 110mm(Diff. 5.0mm)

CANCELLOUS SCREW (32TH) 6.5MM

Catalog no.	011.412.20 to 011.412.110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	20mm to 110mm(Diff. 5.0mm)



CANCELLOUS SCREW (FULL TH) 6.5MM

Catalog no.	011.413.25 to 011.413.110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	25mm to 110mm(Diff. 5.0mm)



CANNULATED SCREW (16TH) 6.5MM

Catalog no.	011.414.20 to 011.414.110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	20mm to 110mm(Diff. 5.0mm)



CANNULATED SCREW (32TH) 6.5MM

Catalog no.	011.415.25 to 011.415.110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	25mm to 110mm(Diff. 5.0mm)





CANNULATED SCREW (FULL THREAD) 6.5MM

Catalog no.	011.417.20 to 011.417.110
Material	Stainless steel 316L,316LVM, Titanium
Size	6.5mm
Length	20mm to 100mm(Diff. 5.0mm)



HARBAT SCREW

Catalog no.	011.416.25 to 011.416.30
Material	Stainless steel 316L,316LVM, Titanium
Size	4.0mm
Length	25mm to 30mm(Diff. 2.0mm)



DHS LAG SCREW

Catalog no.	011.418.50 to 011.418.110
Material	Stainless steel 316L,316LVM, Titanium
Size	12.5mm
Length	50mm to 110mm(Diff. 5.0mm)



DHS LAG SCREW (HEX)

Catalog no.	011.01.418.50 to 011.01.418.110
Material	Stainless steel 316L,316LVM, Titanium
Size	12.5mm
Length	50mm to 110mm(Diff. 5.0mm)

DHS LAG SCREW PEDIATRICS

Catalog no.	011.463.40 to 011.463.105
Material	Stainless steel 316L,316LVM, Titanium
Size	
Length	40mm to 105mm(Diff. 5.0mm)



WASHER FOR CANCELLOUS SCREW

Catalog no.	011.4216
Material	Stainless steel 316L,316LVM, Titanium
Size	4.0mm, 7.0mm, 13mm
Length	Standard





K NAIL FOR FEMUR

Catalog no.	011.545.20 to 011.545.48
Locking Bolt ø	
Material	Stainless steel 316L,316LVM ,Titanium
Nail Diameter	07mm, 08mm, 09mm, 10mm, 11mm, 12mm
Nail Length	20cm to 46cm (Diff. 2.0cm)




V NAIL FOR TIBIA

Catalog no.	011.546.20 to 011.546.40
Locking Bolt ø	
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	07mm, 08mm, 09mm, 10mm, 11mm,
Nail Length	20cm to 40cm (Diff. 2.0cm)



SUPRA CONDYLAR NAIL

Catalog no.	011.548.15 to 011.548.30
Locking Bolt ø	
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	9mm,10mm,11mm,12mm
Nail Length	15cm to30cm(diff.5.0CM)
Proximal & Distal Bolt	
Bolt Length	



RUSH NAIL

Catalog no.	011.534.18 to 011.534.30
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	1.5mm, 2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm, 4.5mm, 5.0mm, 6mm
Nail Length	18cm to 30cm (diff. 2.0cm)

SQUARE NAIL FOR ULNA

Catalog no.	011.547.2.0.18 to 011.547.2.0.30
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm, 4.5mm, 5.0mm
Nail Length	18cm to 30cm (diff. 2.0cm)



SQUARE NAIL FOR RADIUS

Catalog no.	011.547.2.0.18 to 011.547.2.0.30
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm, 4.5mm, 5.0mm
Nail Length	18cm to 30cm (diff. 2.0cm)



HUMERUS INTERLOCKING NAIL

Catalog no.	011.549.18 to 011.549.30
Locking Bolt ø	3.1mm
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	6.0mm, 7.0mm, 8.0mm
Nail Length	18cm to 30cm (diff. 2.0cm)
Proximal Bolt	Locking Bolt 3.1mm
Bolt Length	10mm to 40mm(diff. 2.0mm)
Distal Bolt	Locking Bolt 3.1mm
Bolt Length	10mm to 40mm(diff. 2.0mm)





INTER LOCKING TIBIA NAIL (SOLID)

Catalog no.	011.01.550.26 to 011.01.550.40
Locking Bolt ø	4.9mm
Material	Stainless steel 316L,316LVM ,Titanium
Nail Diameter	8.0mm, 9.0mm,10.0mm,11.0mm
Nail Length	26.0cm to 40.0cm(diff. 2.0cm)
Proximal & Distal Bolt	Locking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)
Proximal & Distal Bolt	Locking Bolt 4.9mm
Bolt Length	65mm to 80mm(diff. 5.0mm)



INTER LOCKING TIBIA NAIL

Catalog no.	011.550.26 to 011.550.40
Locking Bolt ø	4.9mm
Material	Stainless steel 316L,316LVM ,Titanium
Nail Diameter	8.0mm, 9.0mm,10.0mm,11.0mm
Nail Length	26.0cm to 40.0cm(diff. 2.0cm)
Proximal & Distal Bolt	Locking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)
Proximal & Distal Bolt	Locking Bolt 4.9mm
Bolt Length	65mm to 80mm(diff. 5.0mm)



INTER LOCKING FEMUR NAIL

Catalog no.	011.551.34. to 011.551.44
Locking Bolt ø	4.9mm
Material	Stainless steel 316L,316LVM, Titanium
Nail Diameter	8.0mm, 9.0mm,10.0mm,11.0mm,12.0mm
Nail Length	34.0cm to 44.0cm (Diff. 2.0cm)
Proximal & Distal Bolt	Locking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)
Proximal & Distal Bolt	Locking Bolt 4.9mm
Bolt Length	65mm to 80mm(diff. 5.0mm)

TFN NAIL SHORT

Catalog no.	011.5935
Locking Bolt ø	8.0mm, 6.4mm & 4.9mm
Material	Stainless steel 316L, 316LVM Titanium
Nail Diameter	9.0mm, 10.0mm, 11.0mm, 12.0mm
Nail Length	18cm
Degree	130° / 135°
Proximal Bolt	8.0mm
Bolt Length	40mm to 120mm(diff. 5.0mm)
Proximal Bolt	6.4mm
Bolt Length	35mm to 120mm(diff. 5.0mm)
Distal Bolt	Interlocking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)



PFN NAIL SHORT

Catalog no.	011.57
Locking Bolt ø	5.0mm, 3.7mm & 3.0mm
Material	Stainless steel 316L, 316LVM Titanium
Nail Diameter	9.0mm, 10.0mm, 11.0mm, 12.0mm
Nail Length	25cm
Degree	130° / 135°
Proximal Bolt	8.0mm
Bolt Length	40mm to 120mm(diff. 5.0mm)
Proximal Bolt	6.4mm
Bolt Length	35mm to 120mm(diff. 5.0mm)
Distal Bolt	Interlocking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)



PFN NAIL LONG 130° (LEFT & RIGHT)

Catalog no.	011.57.09 to 011.57.12
Locking Bolt ø	3.0mm & 1.5mm
Material	Stainless steel 316L, 316LVM Titanium
Nail Diameter	9.0mm, 10.0mm, 11.0mm, 12.0mm
Nail Length	25cm
Degree	130° / 135°
Proximal Bolt	8.0mm
Bolt Length	40mm to 120mm(diff. 5.0mm)
Proximal Bolt	6.4mm
Bolt Length	35mm to 120mm(diff. 5.0mm)
Distal Bolt	Interlocking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm)
Standard	Left & Right






ELASTIC NAIL (FLEXIBLE NAIL)

Catalog no.	011.5938
Material	Stainless steel 316L,316LVM Titanium
Nail Diameter	2.0mm, 2.5mm, 3.0mm, 3.5mm, 4.0mm, 4.5mm
Nail Length	44.0cm Standard



SIRUS NAIL (LEFT & RIGHT)

Catalog no.	011.559.34 to 011.559.44
Locking Bolt ø	6.5mm & 4.9mm
Material	Stainless steel 316L,316LVM Titanium
Nail Diameter	9.0mm,10.0mm,11.0mm,12.0mm
Nail Length	34cm to 44cm (diff. 2.0cm)
Proximal Bolt	Cannulated / Solid 6.5mm
Bolt Length	30mm to 110m (diff. 5.0mm)
Proximal & Distal Bolt	Locking Bolt 4.9mm
Bolt Length	20mm to 60mm(diff. 2.0mm) 65mm to 80mm(diff. 5.0mm)
Standard	Left & Right



ENDER'S NAIL

Catalog no.	011.553.18 to 011.553.46
Material	Stainless steel 316L,316LVM Titanium
Nail Diameter	3.0mm, 3.5mm, 4.0mm, 4.5mm,
Nail Length	18cm to 46cm (Diff. 1.0cm)

INTER LOCKING BOLT 3.1MM

Catalog no.	011.633.18 to 011.633.50
Thread Diameter	3.1mm
Material	Stainless steel 316L,316LVM ,Titanium
Size	3.1mm
Length	18mm to 50mm (Diff. 2.0mm)



INTER LOCKING BOLT 4.5MM

Catalogue No.	011.01.652.18 to 011.01.652. 90
Thread Diameter	4.5mm
Core Diameter	
Use drills Diameter	3.2mm
For passing hole	
Material	Stainless steel 316L,316LVM, Titanium
Size	4.5mm
Length	18mm to 90mm (Diff. 2.0mm)



INTER LOCKING BOLT 3.9MM

Catalog no.	011.653.18 to 011.653.80
Thread Diameter	3.9mm
Core Diameter	2.7mm
Use drills Diameter	3.2mm
For passing hole	4.0mm
Material	Stainless steel 316L,316LVM, Titanium
Size	3.9mm
Length	18mm to 50mm (Diff. 2.0mm) 55mm to 80mm (Diff. 5.0mm)



INTER LOCKING BOLT 4.9MM

Catalog no.	011.652.18 to 011.652.90
Thread Diameter	4.9mm
Core Diameter	3.8mm
Use drills Diameter	4.0mm
For passing hole	5.1mm
Material	Stainless steel 316L,316LVM, Titanium
Size	4.9mm
Length	18mm to 60mm (Diff. 2.0mm) 65mm to 90mm (Diff. 5.0mm)





PFN BOLT CANNULATED 6.4MM

Catalog no. 011.65.30 to 011.65.120
Thread Diameter
Core Diameter
Use drills Diameter
For passing hole
Material Stainless steel 316L,316LVM, Titanium
Size 6.4mm
Length 30mm to 120mm (Diff. 5.0mm)



PFN BOLT CANCELLOUS 6.4MM

Catalog no. 011.01.65.30 to 011.01.65.120
Thread Diameter
Core Diameter
Use drills Diameter
For passing hole
Material Stainless steel 316L,316LVM, Titanium
Size 6.4mm
Length 30mm to 120mm (Diff. 5.0mm)



PFN BOLT CANNULATED 8.0MM

Catalog no. 011.66.45 to 120
Thread Diameter
Core Diameter
Use drills Diameter
For passing hole
Material Stainless steel 316L,316LVM, Titanium
Size 8.0mm
Length 45mm to 120mm (Diff. 5.0mm)



PFN BOLT CANNULATED 8.0MM (HEX)

Catalog no. 011.01.66.45 to 011.01.66.120
Thread Diameter
Core Diameter
Use drills Diameter
For passing hole
Material Stainless steel 316L,316LVM ,Titanium
Size 8.0mm
Length 45mm to 120mm (Diff. 5.0mm)

THOMPSON HIP PROSTHESIS

Catlog no.	011.754.37 to 011.754.53
Material	Stainless steel 316L, 316LVM, Titanium
Size	37, 39, 41, 43, 45, 47, 49, 51, 53
Length	Standard



AUSTIN MOORE HIP PROSTHESIS

Catlog no.	011.755.37 to 011.755.53
Material	Stainless steel 316L, 316LVM, Titanium
Size	37, 39, 41, 43, 45, 47, 49, 51, 53
Length	Standard



BIPOLAR FIXED HEAD

Catlog no.	011.78.37 to 011.78.53
Material	Stainless steel 316L, 316LVM, Titanium
Size	37, 39, 41, 43, 45, 47, 49, 51, 53
Length	Standard

MODULAR BIPOLAR CUP

Catlog no.	011.7115.37 to 011.7115.53
Material	Stainless steel 316L, 316LVM, Titanium
Size	37, 39, 41, 43, 45, 47, 49, 51, 53
Length	Small, Medium, Large

MODULAR BIPOLAR STEM

Catlog no.	–
Material	Stainless steel 316L, 316LVM, Titanium
Size	Small, Medium, Large
Length	Standard

SHOULDER PROSTHESIS

Catlog no.	011.751.35 to 011.751.51
Material	Stainless steel 316L, 316LVM, Titanium
Size	35, 37, 39, 41, 43, 45, 47, 49, 51
Length	Standard

AUSTIN MOORE HIP PROSTHESIS LONG STEM

Catlog no.	011.01.755.35 to 011.01.755.51
Material	Stainless steel 316L, 316LVM, Titanium
Size	35, 37, 39, 41, 43, 45, 47, 49, 51
Length	Standard



MONO AXIAL SCREW

Catalog no. 011.826.20 to 011.826.50
Material Titanium
Dia 4.5mm, 5.5mm, 6.5mm, 7.5mm
Length 20mm to 50mm(Diff. 5.0mm)



POLY AXIAL SCREW

Catalog no. 011.828.20 to 011.828.50
Material Titanium
Dia 4.5mm, 5.5mm, 6.5mm, 7.5mm
Length 20mm to 50mm(Diff. 5.0mm)



MONO AXIAL SCREW (D/T)

Catalog no. 011.833.20 to 011.833.50
Material Titanium
Dia 4.5mm, 5.5mm, 6.5mm, 7.5mm
Length 20mm to 50mm(Diff. 5.0mm)



POLY AXIAL SCREW (D/T)

Catalog no. 011.834.20 to 011.834.50
Material Titanium
Dia 4.5mm, 5.5mm, 6.5mm, 7.5mm
Length 20mm to 50mm(Diff. 5.0mm)





MONO AXIAL REDUCTION SCREW

Catlog no.	011.827.20 to 011.827.50
Material	Titanium
Size	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)



POLY AXIAL REDUCTION SCREW

Catlog no.	011.829.20 to 011.829.50
Material	Titanium
Size	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)



MONO AXIAL REDUCTION (D/T)

Catlog no.	011.813.20 to 011.813.50
Material	Titanium
Size	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)




POLY AXIAL REDUCTION (D/T)

Catlog no.	011.895.20 to 011.895.50
Material	Titanium
Size	4.5mm, 5.5mm, 6.5mm, 7.5mm
Length	20mm to 50mm(Diff. 5.0mm)


3.5 MM DIA - MONOAXIAL SCREW

Catalog no.	011.830.12 to 011.830.40
Material	Titanium
Dia	3.5mm, 4.2mm, 5.0mm, 6.0mm
Length	12mm to 20mm(Diff. 2.0mm) 25mm to 40mm(Diff. 5.0mm)



3.5MM DIA – POLYAXIAL SCREW

Catalog no.	011.831.12 to 011.831.40
Material	Titanium
Dia	3.5mm, 4.2mm
Length	12mm to 20mm(Diff. 2.0mm) 25mm to 40mm(Diff. 5.0mm)



LAMINAR HOOK ANGLE BLADE

Catalog no.	011.835.01
Material	Titanium
Dia	Standard
Length	Standard



LAMINAR HOOK NARROW BLADE

Catalog no.	011.831.12 to 011.831.40
Material	Titanium
Dia	3.5mm, 4.2mm
Length	12mm to 20mm(Diff. 2.0mm) 25mm to 40mm(Diff. 5.0mm)





PEDICAL HOOK

Catlog no.	011.835.03
Material	Titanium
Size	Standard
Length	Standard



LEMINAR HOOK WIDE BLADE

Catlog no.	011.835.04
Material	Titanium
Size	Standard
Length	Standard



STAPLE

Catlog no.	011.829.01
Material	Titanium
Size	Standard
Length	Standard

CROSS LINK CONNECTORE 5MM X 3MM ROD

Catlog no.	011.838.01
Material	Titanium
Size	Standard
Length	Standard

3.5mm Dia - CONNECTING ROD PRE BEND

Catalog no. 011.823.3.5
Material Titanium
Dia 3.5mm
Length Standard

CONNECTING ROD

Catalog no. 011.832
Material Titanium
Dia 3.0mm, 3.5mm, 5.0mm, 5.5mm
Length 50mm to 450mm


3.5mm Dia - POSTERIOR OCCIPITAL PLATE WITH ROD

Catalog no. 011.854.3.5.100 to 011.854.3.5.300
Material Titanium
Dia 3.5mm
Length 100mm, 200mm, 250mm, 300mm

POSTERIOR OCCIPITAL SCREW

Catalog no. 011.8133.06 to 011.8133.20
Material Titanium
Dia 3.5mm
Length 06mm to 20mm (Diff. 2.0mm)





ANTERIOR CERVICAL CAGE (SQUARE TYPE)

Catlog no. 011.844.4 to 8.14 to 11.844.4 to 8.16
Material Titanium
Size 4.0mm, 5.0mm, 06mm, 7.0mm, 8.0mm
Length 12mmx14mm, 12mmx16mm



EXPANDABLE CAGE

Catlog no. 011.846.20 to 011.846.75
Material Titanium
Size 10mm, 12mm, 14mm, 16mm, 18mm, 20mm
Length 20mm to 75mm (Diff. 5.0mm)



ECO CAGE

Catlog no. 011.845.15 to 011.845.100
Material Titanium
Size 10mm, 12mm, 14mm, 16mm, 18mm, 20mm
Length 15mm to 100mm (Diff. 5.0mm)



PLIF CAGE

Catlog no. 011.847.20 to 011.847.30
Material Titanium
Size 7.0mm, 8.0mm, 9.0mm, 10.0mm, 11.0mm
Length 20mm, 25mm, 30mm

TLIF CAGE (BANANA TYPE)

Catalog no.	011.856.7 to 011.856.11
Material	Titanium
Dia	7.0mm, 8.0mm, 9.0mm, 10.0mm, 11.0mm
Length	4o, 8o, 11o



PLIF CAGE (BULLET TYPE)

Catalog no.	011.857.20 to 011.857.30
Material	Titanium
Dia	8.0mm, 10.0mm, 12.0mm
Length	20mm, 25mm, 30mm



ANTERIOR CERVICAL PLATE (Z TYPE)

Catalog no.	011.842.25 to 011.842.70
Material	Titanium
Dia	Standard
Length	25.0mm,27.5mm,30.0mm,32.5mm, 35.0mm,37.5mm,40.0mm,42.5mm, 45.0mm,47.5mm,50.0mm,55.0mm, 60.0mm,65.0mm,70.0mm



BONE SCREW FOR ACP 4mm TITANIUM

Catalog no.	011.843.10 to 011.843.24
Material	Titanium
Dia	4.0mm
Length	10mm to 24mm (Diff. 2.0mm)



OCCIPITAL LOOPS

Catlog no.	011.855
Material	Titanium
Size	4.0mm
Length	






Get In Touch

 Champion Surgico. Pacca Garrha,
P.O Box 2885, Pakistan.

 +92 52 429 3283

 info@championsurgico.com

 www.championsurgico.com

